





Product Description

Laerdal Stifneck Select Adjustable Collar

Intended use: To assist in cervical stabilization of the adult patient with suspected cervical spine injury. The Laerdal Select is intended to reduce compression of the cervical spine and is to be used as one component of a system spinal immobilisation.

The adjustable extrication collar is incorporated with the following elements:

1. Single piece collar for use up to 24 hours
2. Available in a Blue/Yellow combination with white foam, or in Green with dark foam
3. Four adult sizes in one adjustable collar
4. Flat design for ease of storage
5. Injection moulded High Density Polyethylene  shell with closed cell foam padding (electrostatic expanded, cross linked polyethylene foam )
6. Loop fastener and 5 rivets made of Nylon 
7. Injection-moulded rivets (neck-band), recessed into the foam for comfort
8. Hook fastener injection-moulded for better performance in varying environmental conditions
9. Serrated upper and lower posterior edges facilitate conformity to patient (rear panel)
10. Large rear panel opening for fluid drainage, air flow, cervical spine visualization and palpation
11. Independent locking mechanisms for size adjustment
12. Locking mechanisms provide audible and tactile confirmation of symmetric size selection
13. Sizing windows allow sizing the collar “off patient”, using same procedure as original Stifneck extrication collar
14. Moulded chin piece anchors mandible
15. Radiolucent, MRI/CT scan compatible
16. Tested for Biocompatibility
17. Motion restriction and neutral alignment equivalent to the original Stifneck extrication collar
18. Instructions for use included in each device package
19. In-moulded instructional graphics on rear panel
20. Manufacturer’s name shown on each device
21. The Select collar is CE marked, class I medical device, **93/42/EEC Medical Device Directive**
22. The Select collar is FDA registered and sold in international markets including but not limited to the Americas, EU countries, Australia, Canada and Japan
23. Intended for single use only, “Do not reuse” 
24. Operating temperature: -25°C to 43°C (-13°F to 110°F)
25. Storage temperature: -34°C to 52°C (-30°F to 125°F)
26. Latex free

Physical dimensions:

Size:	Length: (mm / in)	Back Panel Height: (mm / in)	Front Panel Height: (mm / in)	Tracheal Access Hole: (mm ² / in ²)	Rear Palpation Hole: (mm ² / in ²)
<i>Tall</i>	565.15 / 22.25	139.7 / 5.5	143.0 / 5.63	9651.6 / 14.96	5677.4 / 8.8
<i>Regular</i>	565.15 / 22.25	139.7 / 5.5	128.6 / 5.06	7883.9 / 12.22	5677.4 / 8.8
<i>Short</i>	565.15 / 22.25	139.7 / 5.5	114.3 / 4.50	6206.5 / 9.62	5677.4 / 8.8
<i>No-Neck</i>	565.15 / 22.25	139.7 / 5.5	100.1 / 3.94	4645.2 / 7.20	5677.4 / 8.8