



### Paediatric Immobilization

An immobilization approach designed for the special anatomical differences of the paediatric patient

### Laerdal Spinal Immobilization Solutions



#### Laerdal Spineboard Straps

- Features loop-locks and swivel speed-clip options
- 2-piece construction with dual side-release buckles
- 203 cm (80 in.) in length
- Best Strap (SpineSaver Strapping System) also available (not shown)



#### PadPack Alignment Pads

- Provides additional lift for neutral cervical alignment
- Cost-effective, disposable foam pads
- Can be used with HeadBed II or other head immobilization devices

#### Ordering information

##### Spinal System Solution Sets

- 981000 Spinal System Solution Set with extrication straps (includes BackPack, Adult Select SpeedBlocks Starter Pack & 4 loop-locks)
- 981100 Spinal System Solution Set without speed-clip straps (includes BackPack, Adult Select SpeedBlocks Starter Pack & 4 swivel speed-clip straps)

##### Softneck Extrication Collars

- 980010 Softneck Select Collar\*
- 980012 Softneck Select Collar (green)
- 980020 Softneck Red Select Collar
- 980030 Softneck Tall Collar (Red)
- 980040 Softneck Taper Collar (Red)
- 980050 Softneck Short Collar
- 980060 Softneck Flat Neck Collar
- 980070 Softneck Pediatric Collar
- 980080 Softneck Baby Flat Neck Collar
- 980090 Softneck Set in Carry Bag (all sizes)
- 980100 Carry Bag (only)

##### Spineboards/frames

- 982100 BackPack Spineboard (yellow)
- 982110 BackPack Spineboard (green) (single unit)\*\*
- 982120 BackPack Spineboard (green)
- 982130 BackPack Spineboard (green) (single unit)\*\*
- 982140 BackPack Carry Bag
- 701-0000 Laerdal Flat Pad Backboard/Pack (Single Unit)
- 982112 Premium Strap with loop-locks (single unit)
- 982111 Premium Strap with swivel speed-clip ends (single unit)
- 982117 Best Strap (SpineSaver Strapping System) (single unit)

##### Head Immobilization Devices

- 983070 SpeedBlocks Starter Pack\* (Includes 1 Block Set, Universal Base & 1 Strap & Pad Replacement Set)
- 983080 SpeedBlocks Block Set
- 983090 SpeedBlocks Universal Base\*
- 983096 SpeedBlocks Strap and Pad Replacement Set (3 sets)
- 983100 HeadBed II HD (Single Unit)
- 983110 Softneck Select-HeadBed II HD/LumboPack (Single Unit)
- 983120 PadPack Alignment Pads (Single Unit)

##### Training Curriculum/Videos

- 220-00180 Trauma in Spine Care (CD-ROM) Training Curriculum
- 983061 SpeedBlocks Training Video (DVD)
- 983060 SpeedBlocks Training Video (DVD)

\* Also available in case quantities  
\*\* Minimum order 10 boards. Please allow 8-12 weeks delivery lead time

For more information visit [www.laerdal.com](http://www.laerdal.com)



© 2005 Laerdal Medical AG. All Rights Reserved  
98101EN rev B

### Ultimate Immobilization



### Laerdal Spinal Immobilization Solutions

A systematic approach to spinal immobilization

- Extrication Collars
- Spineboards
- Head Immobilization
- Training Curriculum



[www.laerdal.com](http://www.laerdal.com)



### Adult Immobilization

A systematic approach for neutral alignment of the head, neck, torso and pelvis

## Stifneck Extrication Collars

From the original Stifneck Extrication Collars, to the Stifneck Select and new Stifneck Pedi-Select, our Stifneck family offers outstanding motion restriction using the same easy to use, sizing and application techniques. X-ray translucent and MRI/CT compatible, our color family makes perfect sense for use throughout your patient's emergent care.

### Stifneck Select Adjustable Collar

- Four adult collars in one
- Safe, symmetrical sizing
- Oversized tracheal and back panels for continued assessment and access for advanced airway management
- No adhesives or glue-in-molded hooks for Vetro® strap enhances performance in varying environmental conditions
- Flat storage requires minimal space
- Single-piece construction, no assembly



Stifneck Select Adjustable Collar

ALSO available in GREEN

Pedi-adjustable



Stifneck Pedi-Select Adjustable Collar

### Stifneck Pedi-Select Adjustable Collar

- Three paediatric collars in one
- Same easy sizing and application method as its famous parents
- Oversized tracheal and back panels for continued assessment and access for advanced airway management
- No adhesives or glue-in-molded hooks for Vetro® strap enhances performance in varying environmental conditions
- Built-in oxygen mask/air cannula hooks for ease of securing supplemental oxygen devices

### Stifneck Extrication Collars

- The original one-piece Stifneck designed for superior motion restriction
- 6 individual sizes for adults and children
- Distinct color-coding for quick and efficient size identification
- Available individually or in complete sets with carrying bag



Stifneck Extrication Collars



## Total Spinal System Solutions

From the classroom to patient care

For over 20 years, countless healthcare providers have trained with and used Laerdal spinal immobilization products. Responders know that the better they perform, the better the patient outcome. They've learned to trust Laerdal's Spinal System Solutions to improve patient care.

Starting in the classroom, Laerdal's "Issues in Spinal Care" training and education curriculum helps to strengthen the learner's knowledge of spinal anatomy and proper spinal assessment. It challenges the provider to take a systematic approach to spinal immobilization and refines their patient care skills and immobilization techniques.

In the field, Laerdal's spinal solution products help the provider to transition from their classroom experience to patient care. Our family of Stifneck-Extrication Collars, featuring the Stifneck Select and new Stifneck Pedi-Select adjustable collars, offer peer-reviewed performance for outstanding motion restriction. The BaXstrap Spineboard and Laerdal Pedi-Pad are designed for comfort and help facilitate neutral alignment for both adult and paediatric patients. Laerdal's reusable and disposable head immobilization devices, SpeedBlocks and HeadBed II, offer a cost-effective approach to superior head immobilization.

From our individual immobilization products to our complete Spinal System Solutions, we have everything you need to help get your patients back on their feet.

### "Issues in Spinal Care" Training Curriculum Revised Edition!

Presents concepts and techniques of spinal motion restriction in four independent modules: adult, paediatric, aquatics and spinal field role/act.

Order your free "Issues in Spinal Care" CD-Rom today.



BaXstrap

ALSO available in GREEN



Laerdal Pedi-Pad Backboard Pad



SpeedBlocks Head Immobilizer

## Spineboards

Another key factor in proper system spinal immobilization, the Laerdal BaXstrap Backboard and Laerdal Pedi-Pad Backboard Pad assist in maintaining neutral alignment of your patient. Engineered for durability, ease of use and patient comfort, the BaXstrap and Pad-Pad are made to enhance your patient care.

### BaXstrap Spineboard

- Strong lightweight and durable truly the Chairman of the Boards
- Patient weight capacity of 363 kg (800 lbs)
- X-ray translucent for clear view of the spine
- MRI/CT compatible
- Weighs 6.2 kg (13.5 lbs) and measures 40 x 183 x 45 cm (16 x 72 x 2.5 in.)
- One-piece sealed construction with in-molded speed pins
- Ergonomically designed for ease of use and rescuer safety
- Adult and Paediatric strapping points
- Limited Lifetime Warranty

### Laerdal Pedi-Pad Backboard Pad

- Enhances neutral alignment of the paediatric patient
- Anti-Skid foam construction reduces patient movement
- Essential component in system spinal immobilization of paediatric patients

## Head Immobilization Devices

Offering both reusable and disposable head immobilization devices, the SpeedBlocks Head Immobilizer and HeadBed II (not shown) help you to complete the treatment objectives of system spinal immobilization.

### SpeedBlocks Head Immobilizer

- Designed for durability performance and ease of use
- Cost effective, reusable base and blocks with disposable base pad and block foam
- Low profile, open design for clear and easy patient access
- Multi-axis block movement allows for various patient positioning and provides superior motion restriction
- Fits both adult and paediatric patients two years old and up
- X-Ray translucent and MRI/CT compatible
- Fits most spineboards