

Material Information Sheet

Last updated: 19.07.2007

Catalogue number	Product name
980200	Stifneck Pediatric Collar
980100	Stifneck Baby No-Neck Collar

Part/ Group of parts	Material	CAS-no (if available)	Comments
Neckband	High Density Polyethylene		Yellow Color
Chin Support	High Density Polyethylene		Yellow Color
Rivet, Male	Nylon		
Hook Collardot	Nylon		
Hook and Loop	Nylon		
Foam pads	Low-Density Polyethylene	9002-88-4	Biocompatible according to ISO 10993-1 "Safety for use with humans, skin contact <24 h"
Carry Bag Stifneck	Polyester fabric with PVC backing		

This product complies with:
California Proposition 65 (20 04 2007)

This product contains no latex

Author:

Eivind Tønnessen
Technical Product Specialist R&D
Laerdal Medical