

**Material Information Sheet**

Last updated: 18.06.2007

<b>Catalogue number</b>	<b>Product name</b>
980010	Stifneck Select Collar
980020	Stifneck Pedi-Select Collar

Part/ Group of parts	Material	CAS-no (if available)	Comments
Chin support	High-Density Polyethylene		Yellow colour
Neckband	High-Density Polyethylene		Blue colour
Rivet, Male	Nylon		
Loop strap	Nylon		
Foam pads	Low-Density Polyethylene	9002-88-4	Biocompatible according to ISO 10993-1 "Safety for use with humans, skin contact <24 h"

This product complies with:  
California Proposition 65 (20 04 2007)

This product contains no latex or phthalates

Author:

Eivind Tønnessen  
 Technical Product Specialist R&D  
 Laerdal Medical